

CDS NextGeneration™

## **Hematology Controls and Calibrators**

### **CDS 3-PD Hematology Controls**

The CDS 3-PD Hematology Control is designed to document and monitor values obtained from CBC and 3-part differential instruments.

Tri-Pak (CDS 502-004)  
(2 X 3.3 mL LOW)  
(2 X 3.3 mL NORMAL)  
(2 X 3.3 mL HIGH)

→ 110-day closed and 14-day open vial stability allows use for 3 QC months

### **Hematology Calibrator**

The CDS Hematology Calibrator is designed to verify instrument calibration.

1 x 3.0 mL (CDS 502-040)

→ 60-day closed vial stability from assay and 5-day open vial stability

→ Assayed for most 3-part DIFF hematology analyzers

→ Quality Assurance Program available on request

*To ensure optimal performance, refer to the package insert control use.*

### **Product Guarantee**

When used as recommended, performance is guaranteed. If it is confirmed that any CDS product caused a performance problem, CDS will reimburse the laboratory for any reasonable costs incurred as a result.

If it is inferred or threatened that your Warranty or Service Agreement may be altered or cancelled due to the use of CDS reagents, contact CDS at 800-453-3328. It is considered "Unlawful Tying" to require the purchase of one product (reagents) to gain access to a second product (service) under the Sherman, Wilson, and Clayton Acts.



**Clinical Diagnostic Solutions, Inc.** 1800 NW 65th Avenue, Plantation, Florida 33313 USA  
Telephone: 954-791-1773 Fax: 954-791-7118 Toll Free: 800-453-3328 Website: [www.cdsolinc.com](http://www.cdsolinc.com)

*Where customers are Still #1*



202081E R 06.22.07